Package 'VisiumIO'

November 7, 2025

Title Import Visium data from the 10X Space Ranger pipeline

Version 1.7.0

Description The package allows users to readily import spatial data obtained from either the 10X website or from the Space Ranger pipeline. Supported formats include tar.gz, h5, and mtx files. Multiple files can be imported at once with *List type of functions. The package represents data mainly as SpatialExperiment objects.

License Artistic-2.0

Depends R (>= 4.5.0), TENxIO

Imports BiocBaseUtils, BiocGenerics, BiocIO (>= 1.15.1), jsonlite, methods, S4Vectors, SingleCellExperiment, SpatialExperiment, SummarizedExperiment

Suggests arrow, BiocStyle, data.table, knitr, readr, rmarkdown, sf, tinytest

biocViews Software, Infrastructure, DataImport, SingleCell, Spatial

VignetteBuilder knitr

Encoding UTF-8

Roxygen list(markdown = TRUE)

RoxygenNote 7.3.3

BugReports https://github.com/waldronlab/VisiumIO/issues

URL https://github.com/waldronlab/VisiumIO

Collate 'TENxGeoJSON.R' 'TENxSpatialCSV.R' 'TENxSpatialList-class.R' 'TENxSpatialParquet.R' 'TENxVisium-class.R' 'TENxVisiumList-class.R' 'TENxVisiumHD-class.R' 'VisiumIO-package.R' 'utilities.R'

Date 2025-10-14

git_url https://git.bioconductor.org/packages/VisiumIO

git branch devel

git_last_commit 459a588

git_last_commit_date 2025-10-29

2 VisiumIO-package

Maintainer Marcel Ramos <marcel.ramos@sph.cuny.edu>

Contents

Index		2
	TENxVisiumList-class	18
	TENxVisiumHD-class	
	TENxVisium-class	12
	TENxSpatialParquet-class	10
	TENxSpatialList-class	
	TENxSpatialCSV-class	1
	TENxGeoJSON-class	
	st_invert_y	4
	compareBarcodes	3
	VisiumIO-package	2

VisiumIO-package

VisiumIO: Import Visium data from the 10X Space Ranger pipeline

Description

The package allows users to readily import spatial data obtained from either the 10X website or from the Space Ranger pipeline. Supported formats include tar.gz, h5, and mtx files. Multiple files can be imported at once with *List type of functions. The package represents data mainly as SpatialExperiment objects.

Author(s)

Maintainer: Marcel Ramos <marcel.ramos@sph.cuny.edu> (ORCID)

Authors:

- Dario Righelli [contributor]
- Helena Crowell [contributor]

See Also

Useful links:

- https://github.com/waldronlab/VisiumIO
- Report bugs at https://github.com/waldronlab/VisiumIO/issues

compareBarcodes 3

comp	arc	Dar	000	مما
COIIII	$AI \leftarrow$	יהמי	('()()	

Compare barcodes between raw and filtered data

Description

This function compares the barcodes between raw and filtered data **depending** on the order of processing. Typically, the "raw" barcodes are compared to the "filtered" ones. The presence of raw barcodes in the filtered data are marked as TRUE in the resulting data.frame.

Usage

```
compareBarcodes(
  from_resource,
  to_resource,
  spacerangerOut,
  format = c("mtx", "h5"),
  processing = c("raw", "filtered"),
  ...
)
```

Arguments

from_resource character(1) The path to the resource file whose barcodes are used as the basis

of the comparison; typically, the "raw" feature barcodes are used.

to_resource character(1) The path to the resource file whose barcodes are compared to the

from_resource; typically, the "filtered" feature barcodes.

spacerangerOut character(1) A single string specifying the path to the directory where the out-

put of spaceranger count is located; typically (but not necessarily), this is the outs directory. The directory must contain the (processing)_feature_bc_matrix

and spatial sub directories.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

processing character(2) A vector of length 2 that corresponds to the processing type. The

processing types are typically "raw" and "filtered". These are the prefixes of the folder names raw_feature_bc_matrix and filtered_feature_bc_matrix. The order of the vector determines the comparison. By default, processing = c("raw", "filtered"), which means barcodes in the raw data are compared

to the filtered data.

... Additional arguments passed to TENxH5 or TENxFileList.

Value

A data frame with barcodes of the first element in the processing data type as the first column and a logical vector indicating whether the barcodes are found in the second element in processing. For example, if processing is c("raw", "filtered"), then the first column will be the barcodes

st_invert_y

in the raw data and the second column will be a logical vector indicating whether the barcodes are found in the filtered data.

Examples

```
if (interactive()) {
    compareBarcodes(
        from_resource = "V1_Adult_Mouse_Brain_raw_feature_bc_matrix.tar.gz",
        to_resource =
            "V1_Adult_Mouse_Brain_filtered_feature_bc_matrix.tar.gz",
   )
    compareBarcodes(
        from_resource =
            "V1_Adult_Mouse_Brain_raw_feature_bc_matrix.h5",
        to_resource =
            "V1_Adult_Mouse_Brain_filtered_feature_bc_matrix.h5"
   )
    compareBarcodes(spacerangerOut = "~/data/outs", format = "h5")
    compareBarcodes(
        spacerangerOut = "~/data/feature_bc_matrix", format = "mtx"
    compareBarcodes(
        spacerangerOut = "~/data/folder_feature_bc_matrix", format = "mtx"
   )
}
```

st_invert_y

Flip the Y-axis of cell or nucleus segmentations to align with H&E image

Description

This function flips the Y-axis of cell or nucleus segmentations stored in an sf object to align with the H&E image. The Y-axis flipping is necessary because the origin (0,0) in image coordinates is at the top-left corner, while in Cartesian coordinates, the origin is at the bottom-left corner. The function takes into account the image height and scaling factor to accurately flip the Y-coordinates of the segmentations.

Usage

```
st_invert_y(sf, type = c("POINT", "POLYGON"), img_height, scalef)
```

TENxGeoJSON-class 5

Arguments

sf	sf an sf class object read from a .geojson file.
type	character(1) "POINT" for cell centroid, or "POLYGON" for cell segmentation mask. Default is "POINT".
img_height	$numeric (1) \ The \ total \ length \ along \ the \ Y \ axis \ of \ the \ image. \ Obtained \ by \ reading \ in \ hires \ or \ lowres \ .png \ under \ /spatial \ folder \ with \ magick: :image_read().$
scalef	<pre>numeric(1) The scaling factor from a /spatial/scalefactors_json.json</pre>

file

Value

an sf object with Y-axis of the points or polygons flipped

Author(s)

Estella YiXing Dong

Examples

```
geojson_file <- system.file(
    file.path("extdata", "segmented_outputs", "cell_segmentations.geojson"),
    mustWork = TRUE, package = "VisiumIO"
)
geo_data <- sf::st_read(geojson_file, quiet = TRUE)
st_invert_y(
    sf = geo_data, type = "POLYGON", img_height = 3886, scalef = 0.079
)</pre>
```

TENxGeoJSON-class

Import 10X Genomics GeoJSON files

Description

TENxGeoJSON is a class to represent and import GeoJSON files from 10X Genomics. It is a composed class of TENxIO::TENxFile.

Usage

```
TENxGeoJSON(resource)
## S4 method for signature 'TENxGeoJSON,ANY,ANY'
import(con, format, text, ...)
```

6 TENxGeoJSON-class

Arguments

character(1) The path to the file resource con The connection from which data is loaded or to which data is saved. If this is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than writing to a connection. format The format of the output. If missing and con is a file name, the format is derived from the file extension. This argument is unnecessary when con is a derivative of BiocFile. If con is missing, this can be a character vector directly providing the string data text to import.

... Additional inputs to the low level class generator functions

Details

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class in the background.

Value

```
TENxGeoJSON(): An object of class TENxGeoJSON import-method: An sf and data.frame with the GeoJSON data
```

See Also

https://www.10xgenomics.com/support/software/xenium-ranger/3.0/analysis/segmentation-inputs

```
segout_folder <- system.file(
    file.path("extdata", "segmented_outputs"),
    package = "VisiumIO"
)
geojsonres <- file.path(segout_folder, "cell_segmentations.geojson")

TENxGeoJSON(geojsonres)

TENxGeoJSON(geojsonres) |>
    import()
```

TENxSpatialCSV-class Represent and import spatial CSV data from 10X Genomics

Description

TENxSpatialCSV is a class to represent and import spatial CSV files with specific column names. It is a composed class of TENxIO::TENxFile and contains additional slots for the column names and whether the CSV is a list-type of file.

Usage

```
TENxSpatialCSV(resource, colnames = .TISSUE_POS_COLS)
## S4 method for signature 'TENxSpatialCSV,ANY,ANY'
import(con, format, text, ...)
```

Arguments

resource	character(1) The path to the file
colnames	character() A vector specifying the column names of the CSV, defaults to c("barcode", "in_tissue", "array_row", "array_col", "pxl_row_in_fullres", "pxl_col_in_fullres"). Mainly used for the "positions" CSV type of file which does not include column names in the file.
con	The connection from which data is loaded or to which data is saved. If this is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than writing to a connection.
format	The format of the output. If missing and con is a file name, the format is derived from the file extension. This argument is unnecessary when con is a derivative of BiocFile.
text	If con is missing, this can be a character vector directly providing the string data to import.
	Additional inputs to the low level class generator functions

Details

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class in the background. The column names are set to the default values of c("barcode", "in_tissue", "array_row", "array_col", "pxl_row_in_fullres", "pxl_col_in_fullres"). The column names can be changed by specifying the colnames argument in the constructor function.

Set the option "VisiumIO.csvreader" to either "data.table" or "readr" to use the data.table::fread or readr::read_csv functions, respectively. These options are useful when the CSV file is relatively large and the user wants to use faster read-in options. Note that the outputs will still be converted to DataFrame when incorporated to the SpatialExperiment or SingleCellExperiment object.

Value

```
TENxSpatialCSV: An object of class TENxSpatialCSV import-method: A DataFrame object containing the data from the CSV file
```

Slots

```
isList logical(1) A scalar specifying whether the CSV is a list-type of file colnames character() A vector specifying the column names of the CSV variant character(1) A scalar specifying the variant of the CSV file "positions", "cell_boundaries", or "other". The variant is determined by the name of the CSV file within the constructor function. Values include "positions", "cell_boundaries", and "other".
```

compressed logical(1) A scalar specifying whether the CSV is compressed (mainly with a .gz file extension).

Examples

```
sample_dir <- system.file(
    file.path("extdata", "10xVisium", "section1"),
    package = "VisiumIO"
)
spatial_dir <- Filter(
    function(x) endsWith(x, "spatial"), list.dirs(sample_dir)
)
csvresource <- file.path(spatial_dir, "tissue_positions_list.csv")
TENxSpatialCSV(csvresource)
head(import(TENxSpatialCSV(csvresource)), 4)

import(TENxSpatialCSV(csvresource)) |>
    attr("metadata") |>
    lapply(names)
```

TENxSpatialList-class A class to represent and import spatial Visium data

Description

This class is a composed class of TENxFileList, which can contain a list of TENxFile objects, and a TENxSpatialList object. It is meant to handle spatial Visium data from 10X Genomics.

TENxSpatialList-class

Usage

```
TENxSpatialList(
  resources,
  sample_id = "sample01",
  images = c("lowres", "hires", "detected", "aligned", "aligned_fiducials", "cytassist"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions.*",
  bin_size = character(0L),
  ...
)

## S4 method for signature 'TENxSpatialList,ANY,ANY'
import(con, format, text, ...)
```

Arguments

resources	A TENxFileList object or a file path to the tarball containing the matrix / assay data resources.
sample_id	character(1) A single string specifying the sample ID.
images	character() A vector specifying the images to be imported; can be one or multiple of "lowres", "hires", "detected", "aligned".
jsonFile	character(1) A single string specifying the name of the JSON file containing the scale factors.
tissuePattern	character(1) A single string specifying the pattern to match the tissue positions file.
bin_size	character(1) The bin size of the images to import. The default is 008. It corresponds to the directory name square_000um where 000 is the bin value.
	Parameters to pass to the format-specific method.
con	The connection from which data is loaded or to which data is saved. If this is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than writing to a connection.
con	is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

Details

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class.

Value

A SpatialExperiment object

Methods (by generic)

• import(con = TENxSpatialList, format = ANY, text = ANY): Import a TENxSpatialList object

Slots

```
images character() The image name(s) to use with grep and include in the list of files. Can be one of "lowres", "hires", "lowres", "hires", "detected", "aligned", "aligned_fiducials", or "cytassist".
```

scaleJSON character(1) The file name of the scale factors JSON file, defaults to 'scalefactors_json.json'.

tissuePos character(1) The file name of the tissue positions file; typically a .parquet or .csv file.

sampleId character(1) A scalar specifying the sample identifier.

binSize The bin size of the images to import. The default slot value is character(). It typically corresponds to the directory name square_000um where 000 is the bin value.

Examples

```
spatial_dir <- system.file(
    file.path("extdata", "10xVisium", "section1", "outs", "spatial"),
    package = "VisiumIO"
)

TENxSpatialList(resources = spatial_dir, images = "lowres")

TENxSpatialList(resources = spatial_dir, images = "lowres") |>
    metadata() |> lapply(names)
```

TENxSpatialParquet-class

Represent and import spatial Parquet data from 10X Genomics

Description

TENxSpatialParquet is a class to represent and import spatial Parquet files with specific column names. It is a composed class of TENxIO::TENxFile and contains additional slots for the column names and whether the Parquet is a list-type of file.

Usage

```
TENxSpatialParquet(resource, colnames = .TISSUE_POS_COLS)
## S4 method for signature 'TENxSpatialParquet,ANY,ANY'
import(con, format, text, ...)
```

Arguments

resource character(1) The path to the file

colnames character() A vector specifying the column names of the Parquet, defaults to

c("barcode", "in_tissue", "array_row", "array_col", "pxl_row_in_fullres",

"pxl_col_in_fullres").

con The connection from which data is loaded or to which data is saved. If this

is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

text If con is missing, this can be a character vector directly providing the string data

to import.

... Additional inputs to the low level class generator functions

Details

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class in the background. The column names are set to the default values of c("barcode", "in_tissue", "array_row", "array_col", "pxl_row_in_fullres", "pxl_col_in_fullres"). The column names can be changed by specifying the colnames argument in the constructor function.

Value

```
TENxSpatialParquet(): An object of class TENxSpatialParquet import-method: A DataFrame object containing the data from the Parquet file
```

Slots

colnames character() A vector specifying the column names of the Parquet

```
sample_dir <- system.file(
   file.path("extdata", "binned_outputs", "square_002um", "spatial"),
   package = "VisiumIO"</pre>
```

12 TENxVisium-class

```
pspatial_dir <- Filter(
  function(x) endsWith(x, "spatial"), list.dirs(sample_dir)
)
parquetres <- file.path(spatial_dir, "tissue_positions.parquet")
TENxSpatialParquet(parquetres)
import(TENxSpatialParquet(parquetres))

## metadata in attributes
import(TENxSpatialParquet(parquetres)) |>
  attr("metadata") |>
  lapply(names)
```

TENxVisium-class

A class to represent and import a single Visium Sample

Description

This class is a composed class of TENxFileList which can contain a list of TENxFile objects and a TENxSpatialList object. It is meant to handle a single Visium sample from 10X Genomics.

Usage

```
TENxVisium(
  resources,
  spatialResource,
  spacerangerOut,
  sample_id = "sample01",
  processing = c("filtered", "raw"),
  format = c("mtx", "h5"),
  images = c("lowres", "hires", "detected", "aligned", "cytassist"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions.*\\.csv",
  spatialCoordsNames = c("pxl_col_in_fullres", "pxl_row_in_fullres"),
  ...
)

## S4 method for signature 'TENxVisium,ANY,ANY'
import(con, format, text, ...)
```

Arguments

resources

A TENxFileList object or a file path to the tarball containing the matrix / assay data resources.

spatialResource

A TENxSpatialList object or a file path to the tarball containing the spatial data.

TENx Visium-class 13

spacerangerOut character(1) A single string specifying the path to the directory where the output of spaceranger count is located; typically (but not necessarily), this is the outs directory. The directory must contain the (processing)_feature_bc_matrix and spatial sub directories. sample_id character(1) A single string specifying the sample ID.

character(1) A single string indicating the processing folder available e.g., processing "filtered_feature_barcode_matrix" in the spacerangerOut folder. It can be ei-

ther "filtered" or "raw" (default "filtered"). Only used when spacerangerOut is

specified.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

character() A vector specifying the images to be imported; can be one or images

multiple of "lowres", "hires", "detected", "aligned".

jsonFile character(1) A single string specifying the name of the JSON file containing

the scale factors.

character(1) A single string specifying the pattern to match the tissue positissuePattern

tions file.

spatialCoordsNames

character() A vector of strings specifying the names of the columns in the

spatial data containing the spatial coordinates.

In the constructor, additional arguments passed to TENxFileList; otherwise, not

con The connection from which data is loaded or to which data is saved. If this

> is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

If con is missing, this can be a character vector directly providing the string data text

to import.

Details

Typically, the user will not create an object of this class directly but rather use TENxVisiumList constructor function for multiple samples. Note that the images, jsonFile, tissuePattern, and spatialCoordsNames arguments are only considered when the spacerangerOut argument or both the resources and spatialResource arguments are paths to files.

Value

A SpatialExperiment object

Functions

• import(con = TENxVisium, format = ANY, text = ANY): Import Visium data

14 TENx Visium-class

Slots

```
resources A TENxFileList or TENxH5 object containing the Visium data.

spatialList A TENxSpatialList object containing the spatial

coordNames character() A vector specifying the names of the columns in the spatial data containing the spatial coordinates.

sampleId character(1) A scalar specifying the sample identifier.
```

See Also

https://support.10xgenomics.com/spatial-gene-expression/software/pipelines/latest/output/overview

```
outs_dir <- system.file(</pre>
    file.path("extdata", "10xVisium", "section1", "outs"),
    package = "VisiumIO"
)
## using spacerangerOut folder
tv <- TENxVisium(</pre>
    spacerangerOut = outs_dir, processing = "raw", images = "lowres"
)
import(tv)
## with TENxFileList spacerangerOut input
tvfl <- TENxVisium(</pre>
    spacerangerOut = TENxFileList(outs_dir),
    format = "mtx",
    processing = "raw",
    images = "lowres"
)
import(tvfl)
## check metadata of the object
import(tvfl) |>
    metadata() |>
    lapply(names)
## importing h5 format
tvfl <- TENxVisium(</pre>
    spacerangerOut = outs_dir,
    format = "h5",
    processing = "raw",
    images = "lowres"
)
import(tvfl)
```

TENxVisiumHD-class 15

```
rffolder <- file.path(outs_dir, "raw_feature_bc_matrix")
## using resources and spatialResource inputs
tvfl <- TENxVisium(
    resources = rffolder,
    spatialResource = file.path(dirname(rffolder), "spatial"),
    format = "mtx",
    processing = "raw",
    images = "lowres"
)
import(tvfl)</pre>
```

TENxVisiumHD-class

A class to represent and import multiple Visium HD samples

Description

This class contains a SimpleList of TENxVisiumHD objects each corresponding to one sample. The provided spacerangerOut folder should contain a binned_outputs folder where multiple bin_size subfolders are present, e.g., square_002um.

Usage

```
TENxVisiumHD(
  resources,
  spatialResource,
  spacerangerOut,
  segmented_outputs,
  sample_id = "sample01",
  processing = c("filtered", "raw"),
  format = c("mtx", "h5"),
  images = c("lowres", "hires", "detected", "aligned_fiducials"),
 bin_size = c("008", "016", "002"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions\\.parquet",
  spatialCoordsNames = c("pxl_col_in_fullres", "pxl_row_in_fullres"),
)
## S4 method for signature 'TENxVisiumHD, ANY, ANY'
import(con, format, text, ...)
```

Arguments

resources

A TENxFileList object or a file path to the tarball containing the matrix / assay data resources.

16 TENx VisiumHD-class

spatialResource

A TENxSpatialList object or a file path to the tarball containing the spatial data.

spacerangerOut character(1) A single string specifying the path to the directory where the output of spaceranger count is located; typically (but not necessarily), this is the outs directory. The directory must contain the (processing)_feature_bc_matrix and spatial sub directories.

segmented_outputs

character(1) The path to the segmented_outputs directory

sample_id character(1) A single string specifying the sample ID.

processing character(1) A single string indicating the processing folder available e.g.,

> "filtered_feature_barcode_matrix" in the spacerangerOut folder. It can be either "filtered" or "raw" (default "filtered"). Only used when spacerangerOut is

specified.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

images character() A vector specifying the images to be imported; can be one or

multiple of "lowres", "hires", "detected", "aligned".

bin_size character(1) The bin size of the images to import. The default is 008. It

corresponds to the directory name square_000um where 000 is the bin value.

jsonFile character(1) A single string specifying the name of the JSON file containing

the scale factors.

tissuePattern character(1) A single string specifying the pattern to match the tissue posi-

tions file.

spatialCoordsNames

character() A vector of strings specifying the names of the columns in the

spatial data containing the spatial coordinates.

In the constructor, additional arguments passed to TENxFileList; otherwise, not

used

The connection from which data is loaded or to which data is saved. If this con

is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

If con is missing, this can be a character vector directly providing the string data text

to import.

Details

Typically, the user will provide a path to a directory containing the output of the spaceranger count command. The spaceranger count command outputs a folder containing the "raw" or "filtered" ()_feature_bc_matrix.

Value

A SpatialExperiment object

TENxVisiumHD-class 17

Functions

• import(con = TENxVisiumHD, format = ANY, text = ANY): Import Visium HD data from multiple bin sizes

Author(s)

E. Y. Dong, M. Ramos

```
vdir <- system.file(</pre>
    "extdata", package = "VisiumIO", mustWork = TRUE
## with spacerangerOut folder
TENxVisiumHD(spacerangerOut = vdir, bin_size = "002", images = "lowres")
TENxVisiumHD(spacerangerOut = vdir, bin_size = "002", images = "lowres") |>
    import()
## indicate h5 format
TENxVisiumHD(
    spacerangerOut = vdir, bin_size = "002",
    images = "lowres", format = "h5"
)
TENxVisiumHD(
    spacerangerOut = vdir, bin_size = "002",
    images = "lowres", format = "h5"
) |>
    import()
## use resources and spatialResource arguments as file paths
TENxVisiumHD(
    resources = file.path(
        vdir, "binned_outputs", "square_002um",
        "filtered_feature_bc_matrix.h5"
   ),
    spatialResource = file.path(
        vdir, "binned_outputs", "square_002um",
        "spatial"
   bin_size = "002", processing = "filtered",
    images = "lowres", format = "h5"
) |>
    import()
## provide the spatialResource argument as a TENxFileList
TENxVisiumHD(
    resources = file.path(
        vdir, "binned_outputs", "square_002um",
        "filtered_feature_bc_matrix.h5"
```

18 TENxVisiumList-class

```
),
    spatialResource = TENxFileList(
        file.path(
            vdir, "binned_outputs", "square_002um",
            "spatial"
    bin_size = "002", images = "lowres", format = "h5"
) |>
    import()
seg_outs <- system.file(</pre>
    "extdata", "segmented_outputs", package = "VisiumIO", mustWork = TRUE
TENxVisiumHD(
    segmented_outputs = seg_outs,
    format = "h5",
    images = "lowres"
) |>
    import()
```

TENxVisiumList-class A class to represent and import multiple Visium samples

Description

This class contains a SimpleList of TENxVisium objects each corresponding to one sample.

Usage

```
TENxVisiumList(
  sampleFolders,
  sample_ids,
  processing = c("filtered", "raw"),
  images = c("lowres", "hires", "detected", "aligned"),
  format = c("mtx", "h5"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions.*\\.csv",
  spatialCoordsNames = c("pxl_col_in_fullres", "pxl_row_in_fullres"),
  ...
)

## S4 method for signature 'TENxVisiumList,ANY,ANY'
import(con, format, text, ...)
```

Arguments

sampleFolders character() A vector of strings specifying the directories containing the output of the spaceranger count command.

TENxVisiumList-class 19

sample_ids character() A vector of strings specifying the sample IDs. If not provided, the sample IDs will be the names of the sampleFolders. Therefore, the sample_ids must be the same length as sampleFolders. character(1) A single string indicating the processing folder available e.g., processing "filtered_feature_barcode_matrix" in the spacerangerOut folder. It can be either "filtered" or "raw" (default "filtered"). Only used when spacerangerOut is specified. images character() A vector specifying the images to be imported; can be one or multiple of "lowres", "hires", "detected", "aligned". format The format of the output. If missing and con is a file name, the format is derived from the file extension. This argument is unnecessary when con is a derivative of BiocFile. jsonFile character(1) A single string specifying the name of the JSON file containing the scale factors. tissuePattern character(1) A single string specifying the pattern to match the tissue positions file. spatialCoordsNames character() A vector of strings specifying the names of the columns in the spatial data containing the spatial coordinates. In the constructor, additional arguments passed to TENxFileList; otherwise, not used. The connection from which data is loaded or to which data is saved. If this con is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than writing to a connection. text If con is missing, this can be a character vector directly providing the string data to import.

Details

Typically, the user will provide a path to a directory containing the output of the spaceranger count command. The spaceranger count command outputs a folder containing the "raw" or "filtered" ()_feature_bc_matrix.

Value

A SpatialExperiment object

Functions

import(con = TENxVisiumList, format = ANY, text = ANY): Import multiple Visium samples

20 TENx VisiumList-class

See Also

https://support.10xgenomics.com/spatial-gene-expression/software/pipelines/latest/output/overview

```
sample_dirs <- list.dirs(
    system.file(
        file.path("extdata", "10xVisium"),
        package = "VisiumIO"
    ),
    recursive = FALSE, full.names = TRUE
)

tvl <- TENxVisiumList(
    sampleFolders = sample_dirs,
    sample_ids = c("sample01", "sample02"),
    processing = "raw",
    images = "lowres",
    format = "mtx"
)

import(tvl)</pre>
```

Index

* internal	TENxH5, <i>14</i>
VisiumIO-package, 2	TENxIO::TENxFile, 5, 7, 10
.TENxGeoJSON (TENxGeoJSON-class), 5	TENxSpatialCSV, 8
.TENxSpatialCSV (TENxSpatialCSV-class),	TENxSpatialCSV (TENxSpatialCSV-class), 7
7	TENxSpatialCSV-class, 7
.TENxSpatialList	TENxSpatialList, 8, 12, 14, 16
(TENxSpatialList-class), 8	TENxSpatialList
.TENxSpatialParquet	(TENxSpatialList-class), 8
(TENxSpatialParquet-class), 10	TENxSpatialList-class, 8
.TENxVisium (TENxVisium-class), 12	TENxSpatialParquet, <i>11</i>
.TENxVisiumHD (TENxVisiumHD-class), 15	TENxSpatialParquet
.TENxVisiumList (TENxVisiumList-class),	(TENxSpatialParquet-class), 10
18	TENxSpatialParquet-class, 10
	TENxVisium, 18
BiocFile, 3, 6, 7, 9, 11, 13, 16, 19	TENxVisium (TENxVisium-class), 12
company Damaedaa 2	TENxVisium(), 6, 7, 9, 11
compareBarcodes, 3	TENxVisium-class, 12
<pre>import,TENxGeoJSON,ANY,ANY-method</pre>	TENxVisiumHD, <i>15</i>
(TENxGeoJSON-class), 5	TENxVisiumHD (TENxVisiumHD-class), 15
import, TENxSpatialCSV, ANY, ANY-method	TENxVisiumHD-class, 15
(TENxSpatialCSV-class), 7	TENxVisiumList, <i>13</i>
import, TENxSpatialList, ANY, ANY-method	<pre>TENxVisiumList (TENxVisiumList-class),</pre>
(TENxSpatialList-class), 8	18
import, TENxSpatialParquet, ANY, ANY-method	TENxVisiumList-class, 18
(TENxSpatialParquet-class), 10	
import, TENxVisium, ANY, ANY-method	VisiumIO(VisiumIO-package), 2
(TENxVisium-class), 12	VisiumIO-package, 2
import, TENxVisiumHD, ANY, ANY-method	
(TENxVisiumHD-class), 15	
import, TENxVisiumList, ANY, ANY-method	
(TENxVisiumList-class), 18	
(
SpatialExperiment, 13, 16, 19	
st_invert_y,4	
TENxFile, 8, <i>12</i>	
TENxFileList, 8, 9, 12–16, 19	
TENxGeoJSON, 6	
TENxGeoJSON (TENxGeoJSON-class), 5	
TENxGeoJSON-class, 5	