Package 'DNABarcodeCompatibility'

June 12, 2025

```
Type Package
```

Title A Tool for Optimizing Combinations of DNA Barcodes Used in Multiplexed Experiments on Next Generation Sequencing Platforms

Version 1.25.0

Date/Publication 2025-06-12

Maintainer Céline Trébeau < ctrebeau@pasteur.fr>

Description The package allows one to obtain optimised combinations of DNA barcodes to be used for multiplex sequencing. In each barcode combination, barcodes are pooled with respect to Illumina chemistry constraints. Combinations can be filtered to keep those that are robust against substitution and insertion/deletion errors thereby facilitating the demultiplexing step. In addition, the package provides an optimiser function to further favor the selection of barcode combinations with least heterogeneity in barcode usage.

```
License file LICENSE
Encoding UTF-8
LazyData true
RoxygenNote 6.1.1
Imports dplyr, tidyr, numbers, purrr, stringr, stats, utils, methods,
     Rcpp (>= 0.11.2), BH
LinkingTo Rcpp, BH
Depends R (>= 3.6.0)
Suggests knitr, rmarkdown, BiocStyle, testthat
VignetteBuilder knitr
biocViews Preprocessing, Sequencing
URL https://dnabarcodecompatibility.pasteur.fr/
BugReports https://gitlab.pasteur.fr/ida-public/dnabarcodecompatibility/-/issues
NeedsCompilation yes
git_url https://git.bioconductor.org/packages/DNABarcodeCompatibility
git_branch devel
git_last_commit f5db668
git_last_commit_date 2025-04-15
Repository Bioconductor 3.22
```

2 distance

```
Author Céline Trébeau [cre] (ORCID: <a href="https://orcid.org/0000-0001-6795-5379">https://orcid.org/0000-0001-6795-5379</a>),
Jacques Boutet de Monvel [aut] (ORCID: <a href="https://orcid.org/0000-0001-6182-3527">https://orcid.org/0000-0001-6182-3527</a>),
Fabienne Wong Jun Tai [ctb],
Raphaël Etournay [aut] (ORCID: <a href="https://orcid.org/0000-0002-2441-9274">https://orcid.org/0000-0002-2441-9274</a>)
```

Contents

	DNABarcodeCompatibility-package
	distance
	distance_filter
	experiment_design
	file_loading_and_checking
	get_all_combinations
	get_random_combinations
	IlluminaIndexes
	IlluminaIndexesRaw
	optimize_combinations
Index	11
muex	11

DNABarcodeCompatibility-package

DNABarcodeCompatibility: to find optimised sets of compatible barcodes with least heterogeneity in barcode usage for multiplex experiments performed on next generation sequencing platforms.

Description

The DNABarcodeCompatibility package provides six functions to load DNA barcodes, and to generate, filter and optimise sets of barcode combinations for multiplex sequencing experiments. In particular, barcode combinations are selected to be compatible with respect to Illumina chemistry constraints, and can be filtered to keep those that are robust against substitution and insertion/deletion errors thereby facilitating the demultiplexing step. In addition, the package provides an optimiser function to further favor the selection of compatible barcode combinations with least heterogeneity in barcode usage.

distance

Calculate distance between two barcodes.

Description

The function calculates the distance between two barcodes. The user may choose one of several distance metrics ("hamming", "seqlev", "levenshtein", "phaseshift").

Usage

distance(sequence1, sequence2, metric=c("hamming", "seqlev", "levenshtein", "phaseshift"), cost_sub

distance_filter 3

Arguments

sequence1	The first sequence (a string)
sequence2	The second sequence (a string)

metric The distance metric which should be calculated.

cost_sub The cost weight given to a substitution.

cost_indel The cost weight given to insertions and deletions.

Value

The distance between the two sequences.

References

Buschmann, T. and Bystrykh, L. V. (2013) Levenshtein error-correcting barcodes for multiplexed DNA sequencing. BMC bioinformatics, 14(1), 272. Available from http://www.biomedcentral.com/1471-2105/14/272.

Levenshtein, V. I. (1966). Binary codes capable of correcting deletions, insertions and reversals. In Soviet physics doklady (Vol. 10, p. 707).

Hamming, R. W. (1950). Error detecting and error correcting codes. Bell System technical journal, 29(2), 147-160.

Examples

```
distance("AGGT", "TTCC")
distance("AGGT", "TTCC", metric="seqlev")
```

distance_filter

Select barcode combinations with error correction properties

Description

Filters a list of barcode combinations for a given distance metric (hamming or seqlev) and threshold in order to produce a list of barcodes satisfying the distance constraints.

Usage

```
distance_filter(index_df, combinations_m, metric, d)
```

Arguments

index_df A dataframe containing barcodes identifiers, corresponding DNA sequences along

with GC content and presence of homopolymers.

combinations_m A matrix of compatible barcode combinations.

metric The type of distance (hamming or seqlev or phaseshift).

d The minimum value of the distance.

Details

The "hamming" distance is suitable for correcting substitution errors. The "seqlev" distance is suitable for correcting both substitution and insertion/deletion errors.

4 experiment_design

Value

A filtered matrix containing the identifiers of the barcodes satisfying the distance constraints.

References

Buschmann, T. 2015. The Systematic Design and Application of Robust DNA Barcodes.

Buschmann, T. 2017. DNABarcodes: an R package for the systematic construction of DNA sample tags. Bioinformatics 33, 920–922.

See Also

```
get_all_combinations, get_random_combinations
```

Examples

```
barcodes <- DNABarcodeCompatibility::IlluminaIndexes
m <- get_all_combinations(barcodes, 2, 4)
distance_filter(barcodes, m, "hamming", 3)</pre>
```

experiment_design

Find a set of barcode combinations with least heterogeneity in barcode usage for single and dual indexing

Description

This function uses the Shannon Entropy to identify a set of compatible barcode combinations with least heterogeneity in barcode usage, in the context of single and dual indexing. It performs either an exhaustive or a random-greedy search of compatible DNA-barcode combinations depending on the size of the DNA-barcode population and of the number of samples to be multiplexed.

Usage

```
experiment_design(file1, sample_number, mplex_level, platform = 4,
file2 = NULL, export = NULL, metric = NULL, d = 3, thrs_size_comb,
max_iteration, method)
```

Arguments

file1	The input data file that contains 2 columns separated by a space or a tabulation, namely the sequence identifiers and corresponding DNA sequence
sample_number	Number of libraries to be sequenced.
mplex_level	The number at which the barcodes will be multiplexed.
platform	An integer representing the number of channels (1, 2, 4) of the desired Illumina platform: 1 for iSeq; 2 for NextSeq, NovaSeq, MiniSeq; 4 for HiSeq and MiSeq. 0 represents any other platform than Illumina.
file2	The input data file that contains 2 columns separated by a space or a tabulation, namely the sequence identifiers and corresponding DNA sequence; used for dual-indexing, see details below.
export	If not NULL, results are saved in a csv file at the specified path.

experiment_design 5

metric The type of distance (hamming or seqlev).

d The minimum value of the distance.

thrs_size_comb The maximum size of the set of compatible combinations to be used for the

greedy optimization.

max_iteration The maximum number of iterations during the optimizing step.

method The choice of the greedy search: 'greedy exchange' or 'greedy descent'.

Details

By specifying the total number of libraries and the number of libraries to be multiplexed, this function returns an optimal set of DNA-barcode combinations to be used for sequencing.

In the case of **single indexing**, one must only provide one input file containing a list of DNA barcodes (file1 argument). The file2 argument being optional, the program runs the optimisation for single indexing if this argument is left empty. The output shows the sample ID with its respective single barcode.

In the case of **dual indexing**, one must provide two input files containing DNA barcodes as two separate sets of barcodes. The program will detect these two files and automatically switch to the 'dual indexing' mode. The program then runs the optimisation for each barcode set separately. The output shows the sample ID with its respective pair of barcodes.

The inputs of the algorithm are a list of n distinct barcodes, the number N of required libraries, and the multiplex level k; N = ak, where a is the number of lanes of the flow cells to be used for the experiment.

• Step 1:

This step consists of identifying a set of compatible barcode combinations. Given the number of barcodes and the multiplex level, the total number of barcode combinations (compatible or not) reads:

 $\binom{n}{k}$

If this number is not too large, the algorithm will perform an exhaustive search and output all compatible combinations of k barcodes. Otherwise, it will proceed by picking up combinations at random, in order to identify a large enough set of compatible barcode combinations.

• Step 2:

Finds an optimized set of barcode combinations in which barcode redundancy is minimized (see details in optimize_combinations)

Value

A dataframe containing compatible DNA-barcode combinations organized by lanes of the flow cell.

Examples

```
write.table(DNABarcodeCompatibility::IlluminaIndexesRaw,
txtfile <- tempfile(), row.names = FALSE, col.names = FALSE, quote=FALSE)
experiment_design(file1=txtfile, sample_number=18, mplex_level=3,
platform=4)</pre>
```

6 get_all_combinations

```
file_loading_and_checking
```

Loading and checking DNA barcodes.

Description

Loads the file containing DNA barcodes and analyze barcode content.

Usage

```
file_loading_and_checking(file)
```

Arguments

file

The input data file that contains 2 columns separated by a space or a tabulation, namely the sequence identifiers and corresponding DNA sequence.

Details

This function loads the DNA barcodes from the input file and checks barcodes for unicity (identifier and sequence), DNA content, and equal size. It also calculates the fraction of G and C relative to A and T, as referred to as "GC content", and it detects the presence of homopolymers of length >= 3.

Value

A dataframe containing sequence identifiers, nucleotide sequence, GC content, presence of homopolymers.

Examples

```
write.table(DNABarcodeCompatibility::IlluminaIndexesRaw,
txtfile <- tempfile(), row.names = FALSE, col.names = FALSE, quote=FALSE)
file_loading_and_checking(txtfile)</pre>
```

get_all_combinations Get all compatible combinations.

Description

Finds the exhaustive set of compatible barcode combinations.

Usage

```
get_all_combinations(index_df, mplex_level, platform)
```

Arguments

index_df A dataframe containing barcodes identifiers, corresponding DNA sequences along

with GC content and presence of homopolymers.

mplex_level The number at which the barcodes will be multiplexed. Illumina recommends

to not multiplex more than 96 libraries.

platform An integer representing the number of channels (1, 2, 4) of the desired Illumina

platform: 1 for iSeq; 2 for NextSeq, NovaSeq, MiniSeq; 4 for HiSeq and MiSeq.

0 represents any other platform than Illumina.

Details

Be aware that the total number of combinations may become prohibitively large for large barcode sets and large multiplexing numbers.

Value

A matrix containing the identifiers of compatible barcode combinations.

See Also

```
get_random_combinations
```

Examples

```
get_all_combinations(DNABarcodeCompatibility::IlluminaIndexes, 2, 4)
```

```
get_random_combinations
```

Get a large set of compatible combinations.

Description

Finds a randomly generated set of at most 1000 combinations of compatible barcodes.

Usage

```
{\tt get\_random\_combinations(index\_df, mplex\_level, platform)}
```

Arguments

index_df A dataframe containing barcodes identifiers, corresponding DNA sequences along

with GC content and presence of homopolymers.

mplex_level The number at which the barcodes will be multiplexed. Illumina recommends

to not multiplex more than 96 libraries.

Platform An integer representing the number of channels (1, 2, 4) of the desired Illumina

platform: 1 for iSeq; 2 for NextSeq, NovaSeq, MiniSeq; 4 for HiSeq and MiSeq.

0 represents any other platform than Illumina.

8 IlluminaIndexes

Details

This function is suited if the total number of possible combinations is too high for an exhaustive search to be possible in a reasonable amount of time.

Value

A matrix containing the identifiers of compatible barcode combinations.

See Also

```
get_all_combinations
```

Examples

```
get_random_combinations(DNABarcodeCompatibility::IlluminaIndexes, 3, 4)
```

IlluminaIndexes

Barcode dataset from Illumina with features.

Description

48 barcodes from Illumina TruSeq Small RNA kits along with percentage in CG content and presence of homopolymers of size >= 3

Usage

IlluminaIndexes

Format

A data frame with 48 rows and 4 variables:

Id barcode identifier

sequence DNA sequence

GC_content percentage of G and C relative to A and T

homopolymer presence of nucleotide repetitions (>= 3) ...

IlluminaIndexesRaw 9

IlluminaIndexesRaw	Barcode dataset from Illumina.
--------------------	--------------------------------

Description

48 barcodes from Illumina TruSeq Small RNA kits

Usage

IlluminaIndexesRaw

Format

A data frame with 48 rows and 2 variables:

V1 barcode identifier

V2 DNA sequence ...

 $\begin{tabular}{ll} \textbf{optimize_combinations} & \textit{Find a set of barcode combinations with least heterogeneity in barcode} \\ & \textit{usage} \\ \end{tabular}$

Description

This function uses the Shannon Entropy to identify a set of compatible barcode combinations with least heterogeneity in barcode usage.

Usage

```
optimize_combinations(combination_m, nb_lane, index_number,
thrs_size_comb, max_iteration, method)
```

Arguments

combination_m	A matrix of compatible barcode combinations.
nb_lane	The number of lanes to be use for sequencing (i.e. the number of libraries divided by the multiplex level).
index_number	The total number of distinct DNA barcodes in the dataset.
thrs_size_comb	The maximum size of the set of compatible combinations to be used for the greedy optimization.
max_iteration	The maximum number of iterations during the optimizing step.

method The choice of the greedy search: 'greedy_exchange' or 'greedy_descent'.

Details

N/k compatible combinations are then selected using a Shannon entropy maximization approach. It can be shown that the maximum value of the entropy that can be attained for a selection of N barcodes among n, with possible repetitions, reads:

$$S_{max} = -(n-r)\frac{\lfloor N/n \rfloor}{N}\log(\frac{\lfloor N/n \rfloor}{N}) - r\frac{\lceil N/n \rceil}{N}\log(\frac{\lceil N/n \rceil}{N})$$

where r denotes the rest of the division of N by n, while

$$\lfloor N/n \rfloor$$

and

$$\lceil N/n \rceil$$

denote the lower and upper integer parts of N/n, respectively.

Case 1: number of lanes < number of compatible DNA-barcode combinations

This function seeks for compatible DNA-barcode combinations of highest entropy. In brief this function uses a randomized greedy descent algorithm to find an optimized selection. Note that the resulting optimized selection may not be globally optimal. It is actually close to optimal and much improved in terms of non-redundancy of DNA barcodes used, compared to a randomly chosen set of combinations of compatible barcodes.

Case 2: number of lanes >= number of compatible DNA-barcode combinations

In such a case, there are not enough compatible DNA-barcode combinations and redundancy is inevitable.

Value

A matrix containing an optimized set of combinations of compatible barcodes.

See Also

 $\verb"get_all_combinations", \verb"get_random_combinations", \verb"experiment_design" \\$

Examples

```
m \leftarrow get\_random\_combinations(DNABarcodeCompatibility::IlluminaIndexes, 3, 4) optimize_combinations(m, 12, 48)
```

Index

```
* datasets
    IlluminaIndexes, 8
    IlluminaIndexesRaw, 9

distance, 2
distance_filter, 3
DNABarcodeCompatibility-package, 2
experiment_design, 4, 10

file_loading_and_checking, 6
get_all_combinations, 4, 6, 8, 10
get_random_combinations, 4, 7, 7, 10

IlluminaIndexes, 8
IlluminaIndexes, 8
IlluminaIndexesRaw, 9

optimize_combinations, 5, 9
```